


## BSQ Fume Hood



- safer to use anti-corrosive water tap • microprocessor control system, LED display •
- memory function in case of power-failure • porcelain white PP, resistant to acid, alkali and anti-corrosion •
- front window made of thick transparent toughened glass maximize light and visibility inside the fume hood, providing a bright and open working environment •



# BSQ Fume Hood

<b>Specifications (mm)</b>	1200*850*2350   1500*850*2350   1800*850*2350
<b>Top</b>	19/25mm thick epoxy resin worktop colour (black)
<b>External body</b>	1.2 cold-rolling steel with epoxy powder coated internal and movable guide plate made of 6mm thick phenolic resin panel, based on steel body with glass for safety.
<b>Window</b>	6mm thickness tempered safety glass
<b>Lighting equipment</b>	2*20W antiseptis daylight lamp, lamp house is insulated with gas. Lightness on table not lower than 800LUX
<b>Sink</b>	cold water outlet with cup sink, and 250mm diameter PP air outlet
<b>Socket</b>	four 10A, 220v, electric socket
<b>Operation panel</b>	lighting switch, ON/OFF switch for fan, ON/OFF switch for power